

DATA SHEET

Product Name	:	Recombinant Human Thioredoxin1 (rhTRX1)
Synonyms	:	TXN, TRDX, TRX, TRX1, thioredoxin
Description	:	Thioredoxins are proteins that act as antioxidants by facilitating the reduction of other proteins by cysteine thiol-disulfide exchange. Thioredoxins are found in nearly all known organisms and are essential for life in mammals. Thioredoxin is a 12-kD oxidoreductase enzyme containing a dithiol-disulfide active site. It is ubiquitous and found in many organisms from plants and bacteria to mammals. Two cysteines, CXXC motif, are the key to the ability of thioredoxin to reduce other proteins. Thioredoxin proteins also have a characteristic tertiary structure termed the thioredoxin fold. (Ref. From Wikipedia)
NCBI Accession No.	:	NM_003329.3
Amino acid sequence	:	MVKQIESKTAFQEALDAAGDKLVVVDFSATWCGPCKMIKPFFHSLSEKYSNV IFLEVDVDDCQDVASECEVKCMPTFQFFKKGQKVGEFSGANKEKLEATINEL V*
Molecular Mass	:	11.71 kDa (105 aa)
Protein Tags	:	No tagging
Source	:	E. coli.
Cat. No.	:	JW-H039-0500, JW-H039-1000
Storage	:	Should be at \leq -70 °C as undiluted aliquots of handy size. Avoid repeated freezing and thawing.
Cross Reactivity	:	Human, Mouse, Pig, Rabbit, Rat, Sheep, Chicken

Quality Control

Test items	Specifications
Appearance Purity Specificity	 Clear, colorless liquid Greater than 95 % by SDS-PAGE Using Western blot, detection
Concentration	: 0.1 mg/mℓ, Bradford method
Biological Activity	Measured by its ability to catalyze the reduction of insulin. The reaction leads to precipitation, which can be measured by absorbance at 650nm. The specific activity is >10 A650/min/mg
Endotoxin	: Less than 0.5 EU/ μ g as determined by the LAL method
Formulation	: PBS, 0.1mM DTT (pH7.4) without preservative or carrier proteins.
Stability	: Stable for up to 12 months at -70 °C. Stable for a month at 4 °C.
Sterility	Sterilized through a 0.2 μ m membrane filter and packaged aseptically. Culture for 2 weeks, no growth

Manufacturer

: JW CreaGene Inc.

FL2, Jungang induspia V, 138-6, Sangdaewon-dong, Jungwon-gu, Seongnam-si, Gyeonggi-do, South Korea (462-120) Tel: 82-31-737-3310, Fax: 82-31-737-3301